

Lab Energy Challenge



Reducing the total time the freezer door is open by just five minutes a day saves a significant amount of energy and reduces the risk of sample freeze-thaw.

Minimize Door Opening Tracking Sheet

Minimizing freezer door opening will help save energy and protect samples. Fill out the table to see how much your ULT's door is open. Think about how you can minimize.

1. Place this table on ULT for one week
2. Fill out the table by indicating the time and duration the ULT door is open
3. Note the temperature change when the door is closed again

Date Range:		Team Name:				ULT Location:							
Sunday		Monday		Tuesday		Wednesday		Thursday		Friday		Saturday	
Time Open/ duration	Temp when Open and Closed	Time Open/ duration	Temp when Open and Closed	Time Open/ duration	Temp when Open and Closed	Time Open/ duration	Temp when Open and Closed	Time Open/ duration	Temp when Open and Closed	Time Open/ duration	Temp when Open and Closed	Time Open/ duration	Temp when Open and Closed

Green Labs Program

green.labs@ubc.ca | sustain.ubc.ca/greenlabs